



CONTACT US:

Bio3 Implants GmbH
Hanauer Str.1 - 5, 75181,
Pforzheim, Deutschland

+49 723 160 84 104

www.bio3-implants.com

info@bio3-implants.com

Follow us:



BIO3 GUIDE SURGICAL KIT®
for
TEMPLATE SURGERY




German Quality



GSK

Bio3 Guide Surgical Kit (GSK) is a set of drills and instruments required for implant placement during elective surgery, using specialized software and a template. The navigation template is modeled in the software environment based on computed tomography images and 3D scans.

START DRILLS



DG-1.9/8

DG-1.9/10

DG-1.9/11.5

DG-1.9/13

Start Drill D1.9 – Starting drills are used after the First Drill to set the depth and direction for subsequent drilling. Drill D1.9x13 has a total length of 15 mm.








DG-1.9/6

First Drill D1.9 L6 – This is the first drill used to form the initial opening in the bone and set the direction for subsequent drilling.

FINISH DRILLS

Finish Drills are used for the final formation of the bone bed for the implant.



Tissue Punch is a mucotome used to perforate soft tissue through a navigation template. Tissue Punch has a surgical handpiece connection.



G55/TP5.0 (TP5.0)

Mount Driver is used for connecting Moulder Platform with a surgical handpiece.



MD Hex.4

Bio3 GUIDE SURGICAL KIT PRODUCT COMPONENTS

Moulder Platform Short & Moulder Platform Long

Moulders are essential for extracting the implant from the container, placing it into the prepared bone socket, and finally tightening the implant. Moulders are fixed to the implant via screw fixation. They come with adapters for surgical handpieces (Mount Driver) and for a ratchet wrench (Ratchet Adaptor).



Bio3 GUIDE SURGICAL KIT PRODUCT COMPONENTS

Ratchet Adaptor is a wrench for connecting Moulder Platform with a ratchet wrench. GSK includes a short and a long adaptor.



Extractor Ex-M2.5 is used to unfasten implant from Moulder platform in case it is jammed.



Pin for Template is a pin for fixing the navigation template to patient's jaw.



Stop Ring H=2 mm is a set of plastic locking rings 2 mm high, for adjusting the offset. Bio3 GSK includes 4 rings.



Pin Drill D2 is a drill bit for preparing an opening in the bone when fixing the navigation template to patient's jaw.



ADDITIONAL COMPONENTS (not included)

Sleeve for Template (not included in the standard GSK) is a guide bushing for the Pin Drill D2 and Pin for Template that is inserted into the navigation template. Sleeve for Template sets the direction for Pin Drill D2.





Drill Guide DG-D5 and DG-D5-h3.5 (not included in the standard GSK) are guide bushings, which are inserted into the navigation template and set the direction for drills.





TECHNICAL PARAMETERS OF DRILL GUIDES

Drill Guide DG-D5



Drill Guide DG-D5-h3.5



OFFSET SELECTION

Offset is the distance between the neck of the implant and the occlusal surface of the lip of the Drill Guide DG-D5 or Drill Guide DG-D5-h3.5.

Using Offset gives you the ability to adjust the drill position during the procedure for the desired result.


Note that 13 mm drills have an actual length of 15 mm, which can be compensated for using a stop ring or offset.

The maximum implant length that can be placed with this navigation kit is 13 mm and the diameter is 5.0 mm.


Drill Guide DG-D5 or Drill Guide DG-D5-h3.5 have variable Offset values. The table below shows the Offset values according to the implant lengths.

Implant length	Offset value, mm		
L = 6 mm	8.5	10.5	12.5
L = 8 mm	8.5	10.5	12
L = 10 mm	8.5	10	11.5
L = 11.5 mm	8.5	10	12
L = 13 mm	8.5	10.5	-

Installation of Drill Guide and implant:



OFFSET SELECTION



How to determine a rotation marker and ensure the correct orientation of the implant?

The flat side of the hex must be parallel to the flat side of the drill guide (index).

Mounter Platform Short

is used for implant placement, with a standard offset value of 8.5 mm.



Mounter Platform Long

is used for implant placement with an offset value more than 8.5 mm. It has laser markings for offset levels.














Depending on the selected offset value, different drill lengths are used. The table below shows the drill selection depending on the offset value and the implant length.














Implant length	Offset value, mm					
	8.5	10	10.5	11.5	12	12.5
L = 6 mm	Drill L8 (Ring 2 mm)	-	Drill L8	-	-	Drill L10
L = 8 mm	Drill L8	-	Drill L10	-	Drill L11,5	-
L = 10 mm	Drill L10	Drill L11.5	-	Drill L15 (Ring 2 mm)	-	-
L = 11.5 mm	Drill L11.5	Drill L15 (Ring 2 mm)	-	-	Drill L15	-
L = 13 mm	Drill L15 (Ring 2 mm)	-	Drill L15	-	-	-

DRILLING PROTOCOL

Protocol for placement of a **6 mm** implant with offset 8.5 mm

<p>Tissue Punch G55/TP5.0 + Drill Guide DG-D5-h3.5</p> 	<p>Bone Mill BMG 4.8 + Drill Guide DG-D5-h3.5</p> 	<p>First Drill D1.9 L6 + Stop ring SR-H2 mm + Drill Guide DG-D5-h3.5</p> 	<p>Start Drill D1.9 L8 + Stop ring SR-H2 mm + Drill Guide DG-D5-h3.5</p> 	<p>Start Drill D2.8 L8 + Stop ring SR-H2 mm + Drill Guide DG-D5-h3.5</p> 	<p>Finish Drill D3.2 L8 + Stop ring SR-H2 mm + Drill Guide DG-D5-h3.5</p> 	<p>Finish Drill D3.6 L8 + Stop ring SR-H2 mm + Drill Guide DG-D5-h3.5</p> 	<p>Finish Drill D4.0 L8 + Stop ring SR-H2 mm + Drill Guide DG-D5-h3.5</p> 	<p>PT14.2/8 + Mounter Platform MP Wide + Stop ring 2 mm SR-H2 mm + Drill Guide DG-D5-h3.5 + *Ratchet adaptor SAD Hex.4 / **Mount Driver MD Hex.4</p>  <p>*for Ratchet adaptor wrench **for hand piece</p>	<p>Finish Drill D4.8 L8 + Stop ring SR-H2 mm + Drill Guide DG-D5-h3.5</p> 	<p>PT15.0/8 + Mounter Platform MP Wide + Stop ring 2 mm SR-H2 mm + Drill Guide DG-D5-h3.5 + *Ratchet adaptor SAD Hex.4 / **Mount Driver MD Hex.4</p>  <p>*for Ratchet adaptor wrench **for hand piece</p>
---	---	--	--	--	---	---	---	---	---	---

Protocol for placement of a **11.5 mm** implant with offset 8.5 mm

<p>Tissue Punch G55/TP5.0 + Drill Guide DG-D5-h3.5</p> 	<p>Bone Mill BMG 4.8 + Drill Guide DG-D5-h3.5</p> 	<p>First Drill D1.9 L6 + Drill Guide DG-D5-h3.5</p> 	<p>Start Drill D1.9 L11.5 + Drill Guide DG-D5-h3.5</p> 	<p>Start Drill D2.8 L11.5 + Drill Guide DG-D5-h3.5</p> 	<p>Finish Drill D3.2 L11.5 + Drill Guide DG-D5-h3.5</p> 	<p>PT13.3/11.5 + Mounter Platform MP Standard + Drill Guide DG-D5-h3.5 + *Ratchet adaptor SAD Hex.4 / **Mount Driver MD Hex.4</p>  <p>*for Ratchet adaptor wrench **for hand piece</p>	<p>Finish Drill D3.6 L11.5 + Drill Guide DG-D5-h3.5</p> 	<p>PT13.8/11.5 + Mounter Platform MP Standard + Drill Guide DG-D5-h3.5 + *Ratchet adaptor SAD Hex.4 / **Mount Driver MD Hex.4</p>  <p>*for Ratchet adaptor wrench **for hand piece</p>	<p>Finish Drill D4.0 L11.5 + Drill Guide DG-D5-h3.5</p> 	<p>PT14.2/11.5 + Mounter Platform MP Standard + Drill Guide DG-D5-h3.5 + *Ratchet adaptor SAD Hex.4 / **Mount Driver MD Hex.4</p>  <p>*for Ratchet adaptor wrench **for hand piece</p>	<p>Finish Drill D4.8 L11.5 + Drill Guide DG-D5-h3.5</p> 	<p>PT15.0/11.5 + Mounter Platform MP Standard + Drill Guide DG-D5-h3.5 + *Ratchet adaptor SAD Hex.4 / **Mount Driver MD Hex.4</p>  <p>*for Ratchet adaptor wrench **for hand piece</p>
---	---	---	--	--	---	--	---	--	---	--	---	--

DRILLING PROTOCOL

Protocol for placement of a **13 mm** implant with offset 11.5 mm

Tissue Punch **G55/TP5.0**
+ Drill Guide **DG-D5-h3.5**



Bone Mill **BMG 4.8**
+ Drill Guide **DG-D5-h3.5**



First Drill **D1.9 L6**
+ Stop ring **SR-H2 mm**
+ Drill Guide **DG-D5-h3.5**



Start Drill **D1.9 L10**
+ Drill Guide **DG-D5-h3.5**



Start Drill **D1.9 L13**
+ Stop ring **SR-H2 mm**
+ Drill Guide **DG-D5-h3.5**



Finish Drill **D2.8 L10**
+ Drill Guide **DG-D5-h3.5**



Finish Drill **D2.8 L13**
+ Stop ring **SR-H2 mm**
+ Drill Guide **DG-D5-h3.5**



Finish Drill **D3.2 L13**
+ Stop ring **SR-H2 mm**
+ Drill Guide **DG-D5-h3.5**



Using Bio3 GSK you can install the following Bio3 implants:

ADVANCED IMPLANTS

Description	Code	D(mm)	L(mm)	Platform
Implant Advanced D3.3 mm L8 mm	AT13.3/8	3.3	8	Standard
Implant Advanced D3.3 mm L10 mm	AT13.3/10	3.3	10	Standard
Implant Advanced D3.3 mm L11.5 mm	AT13.3/11.5	3.3	11.5	Standard
Implant Advanced D3.3 mm L13mm	AT13.3/13	3.3	13	Standard
Implant Advanced D3.8 mm L8 mm	AT13.8/8	3.8	8	Standard
Implant Advanced D3.8 mm L10 mm	AT13.8/10	3.8	10	Standard
Implant Advanced D3.8 mm L11.5 mm	AT13.8/11.5	3.8	11.5	Standard
Implant Advanced D3.8 mm L13mm	AT13.8/13	3.8	13	Standard
Implant Advanced D4.2 mm L6 mm	AT1B4.2/6	4.2	6	Wide
Implant Advanced D4.2 mm L8 mm	AT1B4.2/8	4.2	8	Wide
Implant Advanced D4.2 mm L10 mm	AT1B4.2/10	4.2	10	Wide
Implant Advanced D4.2 mm L11.5 mm	AT1B4.2/11.5	4.2	11.5	Wide
Implant Advanced D4.2 mm L13mm	AT1B4.2/13	4.2	13	Wide
Implant Advanced D5 mm L6 mm	AT1B5/6	5	6	Wide
Implant Advanced D5 mm L8 mm	AT1B5/8	5	8	Wide
Implant Advanced D5 mm L10 mm	AT1B5/10	5	10	Wide
Implant Advanced D5 mm L11.5 mm	AT1B5/11.5	5	11.5	Wide
Implant Advanced D5 mm L13mm	AT1B5/13	5	13	Wide

PROGRESSIVE IMPLANTS

Description
Implant Progressive D3.3 mm L8 mm
Implant Progressive D3.3 mm L10 mm
Implant Progressive D3.3 mm L11.5 mm
Implant Progressive D3.3 mm L13mm
Implant Progressive D3.8 mm L8 mm
Implant Progressive D3.8 mm L10 mm
Implant Progressive D3.8 mm L11.5 mm
Implant Progressive D3.8 mm L13mm
Implant Progressive D4.2 mm L6 mm
Implant Progressive D4.2 mm L8 mm
Implant Progressive D4.2 mm L10 mm
Implant Progressive D4.2 mm L11.5 mm
Implant Progressive D4.2 mm L13mm
Implant Progressive D5 mm L6 mm
Implant Progressive D5 mm L8 mm
Implant Progressive D5 mm L10 mm
Implant Progressive D5 mm L11.5 mm
Implant Progressive D5 mm L13mm

Code	D(mm)	L(mm)	Platform
PT13.3/8	3.3	8	Standard
PT13.3/10	3.3	10	Standard
PT13.3/11.5	3.3	11.5	Standard
PT13.3/13	3.3	13	Standard
PT13.8/8	3.8	8	Standard
PT13.8/10	3.8	10	Standard
PT13.8/11.5	3.8	11.5	Standard
PT13.8/13	3.8	13	Standard
PT1B4.2/6	4.2	6	Wide
PT1B4.2/8	4.2	8	Wide
PT1B4.2/10	4.2	10	Wide
PT1B4.2/11.5	4.2	11.5	Wide
PT1B4.2/13	4.2	13	Wide
PT1B5/6	5	6	Wide
PT1B5/8	5	8	Wide
PT1B5/10	5	10	Wide
PT1B5/11.5	5	11.5	Wide
PT1B5/13	5	13	Wide

FIX IMPLANTS

Description	Code	D(mm)	L(mm)	Platform
Implant FIX D3.5 mm L8 mm	FIX3.5/8	3.5	8	Standard
Implant FIX D3.5 mm L10 mm	FIX3.5/10	3.5	10	Standard
Implant FIX D3.5 mm L11.5 mm	FIX3.5/11.5	3.5	11.5	Standard
Implant FIX D3.5 mm L13mm	FIX3.5/13	3.5	13	Standard
Implant FIX D4 mm L6 mm	FIX4.0/6	4	6	Standard
Implant FIX D4 mm L8 mm	FIX4.0/8	4	8	Standard
Implant FIX D4 mm L10 mm	FIX4.0/10	4	10	Standard
Implant FIX D4 mm L11.5 mm	FIX4.0/11.5	4	11.5	Standard
Implant FIX D4.5 mm L6 mm	FIX4.5/6	4.5	6	Standard
Implant FIX D4.5 mm L8 mm	FIX4.5/8	4.5	8	Standard
Implant FIX D4.5 mm L10 mm	FIX4.5/10	4.5	10	Standard
Implant FIX D4.5 mm L11.5 mm	FIX4.5/11.5	4.5	11.5	Standard
Implant FIX D4.5 mm L13mm	FIX4.5/13	4.5	13	Standard
Implant FIX D5 mm L6 mm	FIX5.0/6	5	6	Standard
Implant FIX D5 mm L8 mm	FIX5.0/8	5	8	Standard
Implant FIX D5 mm L10 mm	FIX5.0/10	5	10	Standard
Implant FIX D5 mm L11.5 mm	FIX5.0/11.5	5	11.5	Standard
Implant FIX D5 mm L13mm	FIX5.0/13	5	13	Standard